

Shodex RSpak Polymer Reversed-Phase Columns

Shodex RSpak Columns are packed with porous polymeric particles that remain stable in a pH range of 2–12. Similar to conventional polymer-based materials, the DS-613 sorbent works well with samples that are more hydrophobic than hydrophilic, and which, consequently, require relatively high concentrations of organic modifiers. DE-613 columns, with a polymethacrylate packing, are more hydrophilic than hydrophobic, and work well with mobile phases containing relatively high concentrations of water. The least hydrophobic sorbent is used for the DE-613 columns.

For weakly cationic species, the DC-613 column is a cation exchanger with unique selectivity (mixed-mode, ion-exchange, and reversed-phase partition chromatography).

Ordering Information

Shodex RSpak D Series Columns

Description	Base-polymer	Functional Group	Dimension	P/N
DS-613	Polystyrene	None	6 × 150 mm	WAT034220
DE-613	Polymethacrylate	None	6 × 150 mm	WAT034221
DC-613	Polystyrene	Sulfonated	6 × 150 mm	WAT034223
DS-G Pre-column	—	—	4.6 × 10 mm	WAT034224
DE-G Pre-column	—	—	4.6 × 10 mm	WAT034225
DC-G Pre-column	—	—	4.6 × 10 mm	WAT034227

Application-Specific Columns

SUGAR AND CARBOHYDRATE ANALYSIS

High-Performance Carbohydrate Analysis Cartridge Column, p/n: WAT044355

Waters High-Performance Carbohydrate Cartridge Column, with reusable end-fittings, is packed with a 4 µm, spherical silica. This column was developed to separate five monosaccharides and disaccharides with baseline resolution in less than 12 minutes. The 4.6 mm I.D. × 250 mm High-Performance Carbohydrate Cartridge Column offers optimal speed, resolution, and longevity. The pre-packed, disposable cartridge column requires reusable end-fittings, which are available separately.

Carbohydrate Analysis Column, p/n: WAT084038

The Carbohydrate Analysis Column uses a covalently bonded amino packing on a silica substrate. It is best suited for low-molecular-weight sugars such as mono-, di-, and tri-saccharides.

Sugar Pak I Column, p/n: WAT085188

The Sugar Pak I Column separates monosaccharides and sugar alcohols via a strong cation-exchange mechanism. The resin is based on a sulfonated styrene-divinylbenzene polymer that provides pH stability by means of a calcium counter ion.

Waters offers a range of columns for the analysis of sugars, carbohydrates, organic acids, and alcohols. Refer to the following tables for ordering information.

Typical Applications for Sugar and Carbohydrate Columns						
Cartridge/Column	Carbohydrate Analysis Column	SAM™ I Reagent with Silica Cartridge	Sugar-Pak™ I, SC-1011, SP-0810	SH-1011, IC-Pak™ Ion-Exclusion Fast Fruit Juice	Dextro-Pak™	KS-800 series
Mode	Partition	Partition	Ion exchange/size exclusion	Ion exchange/size exclusion	Reversed phase	Size exclusion
Eluent	65–85% acetonitrile/water ambient to 70 °C	70–80% acetonitrile/water 0.1% SAM I ambient	Water 75–95 °C	0.01 N phosphoric acid 50–60 °C	Water ambient	—
Application	Mono-, di- and tri-saccharides up to DP 8 sugars and sugar alcohols	Mono-, di- and tri-saccharides	Mono-, di-, oligosaccharides and sugar alcohols	Sugar acids, sugar alcohols, organic acids	Hydrolysed syrups, derivatized sugars	Mono- through oligosaccharides such as syrups
Elution Order	Smallest elute first	Smallest elute first	Largest elute first	Largest and most acidic elute first	Smallest elute first	Largest elute first

Guide to Shodex Sugar Columns

S	C	18	2	1
Type of Column	Cation	% Cross Linkage	Pore Size	0 - Gel Type
S = sugar	H = H ⁺	—	1 = 20Å	1 - Semimacropore gel
	C = Ca ²⁺	—	2 = 50Å	2 - Permanent pore gel
	P = Pb ²⁺	—	3 = 100Å	
	Z = Zn ²⁺	—	4 = 500Å	
	—	—	5 = 1000Å	
Example:				
S	C	10	1	1
Sugar column	Ca ²⁺	10% cross linkage	20Å	Semimacropore gel

Ordering Information

SAM I Reagent Column

Description	Dimension	Qty.	P/N
SAM I Reagent	7.8 × 300 mm	1/pk	WAT010873

Columns for Alcohols and Carbohydrates

Description	Dimension	Particle Size	Qty.	P/N
Carbohydrate Analysis Column	—	—	1/pk	WAT084038
Dextro-Pak Cartridge Column	8.0 × 100 mm	—	1/pk	WAT085650
High-Performance Carbohydrate Sentry Guard Column	3.9 × 20 mm	4 µm	2/pk	WAT046895 ¹
SC-1011 Column	8.0 × 300 mm	—	1/pk	WAT034238
SC-1011P Pre-column	6.0 × 50 mm	—	1/pk	WAT034244
SH-1011	8.0 × 300 mm	—	1/pk	WAT034236
SH-1011P Pre-column	6.0 × 50 mm	—	1/pk	WAT034243
SP-0810 Column	8.0 × 300 mm	—	1/pk	WAT036954
SP-0810P Pre-column	6.0 × 50 mm	—	1/pk	WAT034245
Sugar-Pak 1 Column	6.5 × 300 mm	—	1/pk	WAT085188
Sugar-Pak 1 Guard-Pak Inserts	—	—	10/pk	WAT015209 ²
Shodex KS-801	—	7 µm	1/pk	WAT034276

¹ Requires Sentry Guard Holder, p/n: WAT046905.

² Requires Guard-Pak Holder, p/n: WAT088141.

High-Performance Carbohydrate Analysis Cartridge Column

Description	Dimension	P/N
High-Performance Carbohydrate Cartridge Column (requires end-fittings)	4.6 × 250 mm	WAT044355
Sentry Integrated Guard Holder (for Waters cartridge columns)	—	WAT046905

FERMENTATION ANALYSIS, ORGANIC ACIDS, ALCOHOLS, AND CARBOHYDRATES

The ion-exclusion mode is ideally suited for the separation of monosaccharides, organic acids, or sugar acids. The column packings are sulfonated styrene divinylbenzene resins in the hydrogen form (IC-Pak Ion-Exclusion or SH-1011), and the mobile phase is a dilute acid such as 0.01 N phosphoric acid using column temperatures of 50–60 °C.

In this mode, the Fast Juice column can effectively separate glycerol, acetic acid, and ethanol in grape or other fruit juice. The column can also analyze the degree of microbial defect, the extent of natural fermentation in grapes, and the amount of sulfite in various foods and beverages. The IC-Pak Ion-exclusion Column can separate a wide range of organic acids while the Shodex SH Column separates acids as well as larger carbohydrates.

The analysis of alcohols and organic acids is important, for they typically help determine the flavor characteristics of beverages such as wine, beer, and some distilled spirits. The presence of alcohols in fruit juices can indicate product deterioration. The Shodex KC-811 Column, which provides ion-exchange and reversed-phase chromatography modes, is packed with a sulfonated, rigid, styrene-divinylbenzene copolymer. With high efficiency, this packing separates low-molecular-weight organic acids and water-soluble organics such as alcohols, aldehydes, and nitriles. The column provides ion-exclusion and reversed-phase mode of chromatography. Typical mobile phases, run at 1 mL/min at 45–80 °C, are composed of aqueous solutions containing 1% phosphoric acid, acetic acid, or perchloric acid.

Shodex KC-811 Column Retention Chart for Organic Acids

Sample	Retention Time	Sample	Retention Time
Oxalic Acid	5.20	β- Hydroxy-propionic Acid	8.60
Maleic Acid	5.80	D-Glucuronic Acid	8.65
α-Ketoglutaric Acid	5.90	Fumaric Acid	8.95
Citric Acid	6.20	Formic Acid	9.20
Tartaric Acid	6.55	Acetic Acid	9.80
Pyruvic Acid	6.65	Adipic Acid	9.80
trans-Aconitic Acid	6.95	Levulinic Acid	10.00
Glyoxylic Acid	7.00	Mesaconic Acid	10.40
Malic Acid	7.05	Pyroglutamic Acid	10.70
Malonic Acid	7.07	Propionic Acid	11.25
Citraconic Acid	7.20	Acrylic Acid	11.60
Succinic Acid	8.00	Pivalic Acid	14.05
Glycolic Acid	8.40	Methacrylic Acid	14.10
Itaconic	8.50	trans-Crotonic Acid	15.65
Lactic Acid	8.60		

Eluent: Water with 0.1% phosphoric acid, Temperature: 60 °C, Flow rate: 1 mL/min.

Ordering Information

Columns for Fermentation Analysis, Organic Acids, Alcohols, and Carbohydrates

Description	Dimension	Qty.	P/N
Fast Fruit Juice Analysis	8.0 × 100 mm	1/pk	WAT010639
Fast Fruit Juice Guard-Pak Inserts	—	10/pk	WAT015207 ¹
IC-Pak Ion-Exclusion	7.8 × 300 mm	1/pk	WAT010290
SC-1011 Column	8.0 × 300 mm	1/pk	WAT034238
SC-1011P Pre-column	6.0 × 50 mm	1/pk	WAT034244
KC-811	8.0 × 300 mm	1/pk	WAT034298
KC-811 Pre-column	6.0 × 50 mm	1/pk	WAT035501

¹Requires 7.8 × 10 mm Cartridge Holder, p/n: 186000708.

FREE FATTY ACID ANALYSIS

The Waters Free Fatty Acid HP Column uses a phenyl-bonded packing and a simple isocratic elution method to separate free fatty acids on the basis of carbon-chain length and degree of saturation. The short column dimension (3.9 × 150 mm) significantly reduces analysis time and increases sensitivity.

Column performance is based on:


- Straight chain saturated acids, which elute in order of increasing carbon number
- Unsaturated acids which elute before the analogous saturated compound
- Carbon number and chain configuration: the greater the unsaturation, the earlier the elution

Ordering Information

Free Fatty Acid HP Column

Description	Dimension	Particle Size	Qty.	P/N
Free Fatty Acid HP	3.9 × 150 mm	4 µm	1/pk	WAT011690

Fatty Acid Standards



CARBAMATE ANALYSIS KITS



Waters Carbamate Analysis Kits for environmental and food testing include the Waters Carbamate Column, Oasis HLB Cartridges, vials, and reference standards. When used, in part, with regulated methods, these proven kits simplify your analysis while increasing your confidence in the result.

Ordering Information

Carbamate Analysis Kits

Description	P/N
Carbamate Analysis Kit for Environmental Testing	176001740
Carbamate Analysis Kit for Food Testing	186004719

Carbamate Analysis Column for Pesticides

Description	Dimension	Qty.	P/N
Carbamate Analysis	3.9 × 150 mm	1/pk	WAT035577